


P A E D I A T R Y

OVER 25 YEARS OF KNOW  
HOW DEDICATED TO  
PAEDIATRIC RADIOLOGY

For many years, Balt has offered solutions to the specialists in paediatric interventional radiology, through a wide range of products for angiography, angioplasty, neuroradiology and visceral radiology.

Faithful to its policy, Balt has been working in conjunction with leading physicians for over 25 years, developing devices that

fulfil the most specific requirements and remains responsive to your needs to allow you to benefit from its know how.



DC001 - 2003/09



10, RUE CROIX VIGNERON - 95160 MONTMORENCY - FRANCE

TÉL : (33) (0)1 39 89 46 41 - FAX : (33) (0)1 34 17 03 46

www.balt.fr

P A E D I A T R Y



IVA  
Angiography  
Baltacci  
CBVP  
Nephrostomy



Paediatric  
range



3F  
4F


# IVA

Catheter valve introducer

- The sheath hub is transparent therefore allowing for better control of blood reflux and air bubbles
- The sheath is perfectly adjusted on the dilatator for a non traumatic use
- The improved ergonomics of the introducer facilitates the manipulation
- Sizes 3F and 4F are supplied with a guide of 0.018" (0.46mm)
- Colour coded : 3F white  
4F red



Reference:  
IVA3F  
IVA4F



10, RUE CROIX VIGNERON  
95160 MONTMORENCY - FRANCE  
TÉL : (33) (0)1 39 89 46 41  
FAX : (33) (0)1 34 17 03 46  
www.balt.fr

3F  
4F

# Angiography

Angiography

- Available in three different lengths: 50cm, 70cm or 100cm
- All shapes are available on request, please visit our web site for the complete range
- 3F (guide wire: .021" - Flow rate: 7 ml/s)
- 4F (guide wire: .035" - Flow rate: 10 ml/s)



Examples of references for 3F catheters 70 cm

Reference	Curve
C3F70D	Straight*
C3F70C1	Cobra 1
C3F70JD15	Judkins Right 1,5
C3F70JG15	Judkins Left 1,5
C3F70HE	Head Hunter
C3F70P	Pigtail*
C3F70CE10	Cerebral 10


\*with side holes

1.2FP  
1.5FP  
1.8FP

# Baltacci

Intravascular infusion

- 100cm long, with ORX transition from rigid to supple part for a use without guiding-catheter
- Progressive suppleness and diameter
- DMSO compatible



The maximum use pressure of the BALTACCI is 100 psi (7 bars). To not exceed this pressure (100 psi or 7 bars), it is recommended to use a syringe 2,5 ml minimum.


# Cristal Balloon

Paediatric valvuplasty

CBVP

- 70cm long, small shaft
- Semi-compliant balloon

Reference	Balloon ØxL(mm)	Catheter Shaft	working pressure	Guide wire	Balt IVA
CBVP8x20/70	8x20	5.5F	8 bar	.035"	6F
CBVP10x20/70	10x20	5.5F	8 bar	.035"	6F
CBVP12x20/70	12x20	6F	8 bar	.035"	7F
CBVP12x40/70	12x40	6F	8 bar	.035"	7F
CBVP15x25/70	15x25	7F	6 bar	.035"	8F
CBVP15x40/70	15x40	7F	6 bar	.035"	8F



Many other sizes available

5F

# Nephro

Nephrostomy

The nephrostomy set is composed of:

- a trocar needle and a black Teflon X-Ray 5F x 29cm catheter
- a flexible guide wire .035" x 60cm with a small curve
- a fixed core guide wire .038" x 60cm with a straight flexible tip
- a trocar needle and a Pigtail 5F x 25cm catheter
- a flexible 25cm connecting tube for drainage bag
- a one way stopcock

Reference : SNEPHR05F

